Temperature of reaction chamber	Tradalatin 9.11 G
	Operating Temperature 25°C-75°C, Adjustable
	Specific to ±0.5°C
Heating rate	≥4°C/min
Cooling rate	≥2°C/min
Pump rate	70±20ml/min
Leak Tightness	Well-fixed separator, Completely separated, No leakage
Mode	Manual mode / Program mode
	Program mode: Parameters set can be saved and read
Dimension (L x W x H)	361 mm x 245 mm x 167 mm
Maximum Noise produced	≤65dB
Power	AC~220V/50Hz or AC~110V/50Hz
	Input Power:≤350VA

Resolution: 0.1°C

Hybribio Limited

25/F, Wing Hang Finance Centre, 60 Gloucester Road, Wanchai, Hong Kong

www.hybribio.com