



# VES<sup>TM</sup>MATIC EASY

Automatic instrument for the determination of the **ESR** on 10 samples simultaneously or in random access mode.



**NIECCE  
NIECCE**

INNOVATIVE CLINICAL DIAGNOSTIC SYSTEMS



## VES-MATIC EASY:

*Complete with printer.  
USB connector.  
Ref# 10376*



# technical characteristics

- Up to 10 samples using Ves-tec or Vacu-tec tubes
- Random access mode
- Printed results
- Results in 20 minutes
- External barcode reader connector
- LCD screen of 16 characters per line
- USB port for host connection
- Dimension in centimeters (15 x 15 x 22) (l x w x h)
- External power supply IN100/240 VAC OUT 9VDC
- External barcode reader (optional)

DIESSE INC.

1690 W 38 Place, Unit B1 Hialeah, FL 33012 USA  
Phone: (305) 827-5761 | 1-877-DIESSE-3 | Fax: (305) 827-5762

<http://www.diesse-usa.com> e-mail: [salesoffice@diesse.us](mailto:salesoffice@diesse.us)